



Be Right™



Dissolved Oxygen Reagent Set, Digital Titrator, 300 mL

Product #: 2272200
SGD Price: Contact Hach
Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Dissolved Oxygen by Digital Titrator titration, using the Winkler Azide Modification method. Hach Method 8215. Range: 1 - 10 mg/L Set Includes: Alkaline Iodide-Azide Powder Pillows, Manganous Sulfate Powder Pillows, Sodium Thiosulfate Digital Titrator Cartridge (0.200 N), Starch Indicator Solution and Sulfamic Acid Powder Pillows for approximately 50 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method: 8215
Method Name: Winkler
Number of tests: 50
Platform : Digital Titrator
Range: 1 - 10 mg/L O₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)