



Be Right™



Sodium Thiosulfate Digital Titrator Cartridge, Stabilized, 0.0250 N

Product #: 2409301
SGD Price: Contact Hach
Available

For dissolved oxygen determination by Winkler titration. Digits x 0.1 = mg/L DO in 25 mL. Stabilized sodium thiosulfate formulation

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method Name: Winkler
Platform : Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)