



Be Right™



Ammonia Cyanurate Reagent Powder Pillows, 10 mL, pk/100

Product #:

2653199

SGD Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Ammonia Nitrogen by the Salicylate method. Hach Method 8155.

Range: 0.01 - 0.50 mg/L NH₃-N. Sample Size: 10 mL.

Pack of 100 Powder Pillows.

Specifications

Method Name:	Salicylate
Number of tests:	100
Parameter:	Ammonia
Platform :	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)