



Be Right™



## BOD Standard Solution, 300 mg/L, pk/16 - 10-mL Voluette® Ampules

Product #: 1486510  
SGD Price: Contact Hach  
Available

For accuracy checks in Biochemical Oxygen Demand (BOD) determinations.

300 mg/L Glucose plus Glutamic Acid.

Pack of 16 10-mL ampules.

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Specifications

Application:	Calibration Standard
Concentration:	300 mg/L
Description :	BOD Standard (Glucose plus Glutamic Acid)
Formula:	Glucose & Glutamic Acid (GGA)
Package Type:	Ampuls
Parameter:	BOD
Platform :	Solution
Quantity:	16 pcs
	300 mg/L BOD