



Be Right™



BOD Seed Inoculum, Polyseed, 50 capsules/pk

Product #: 2918700
SGD Price: Contact Hach
Available

BOD Seed Inoculum for use in Biochemical Oxygen Demand (BOD5) test procedures. A blend of specialized microbial cultures in an easy-to-use capsule. Designed to provide a uniform standard for the degradation of both industrial and municipal waste in BOD5 analyses. Package of 50 capsules. International Laboratory Supply (InterLab), LTD

Specifications

Application:	Dilution water preparation
Description :	Seed Source
Formula:	Seed Source
Method:	8043
Method Name:	Dilution
Package Type:	Capsules
Parameter:	Biochemical Oxygen Demand
Quantity:	50 pcs

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)