



Be Right™



## Dissolved Oxygen Reagent Set, Buret, 300 mL sample

**Product #:** 2351400  
**SGD Price:** Contact Hach  
**Call for ship date**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Dissolved Oxygen by the Winkler Method (Azide Modification) and buret titration. Range: up to 10 mg/L as DO. Sample Size: 300-mL. Contains reagents for 100 tests.

---

### Specifications

EPA compliant:	Yes
Method:	8229
Number of tests:	100
Range:	1 - > 10 mg/L O <sub>2</sub>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)