



Be Right™



Filling Solution, Internal, Polarographic DO, 30 mL

Product #: 2759123
SGD Price: Contact Hach

Filling solution for Polarographic Dissolved Oxygen probe, Potassium Chloride/Glycerine, 30 mL bottle

Specifications

Description : 5130, 5197000
Formula: KCl & Glycerine
Manual Languages: Okay. 16.10.2012 Bauer
Quantity: 30 mL
Solution Type: Refill

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)