



Be Right™



## Detergents Reagent Set (No Air Shipments)

Product #:

2446800

 **Hazardous**

SGD Price:

Contact Hach

No availability shown

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Detergents (Surfactants) by the Crystal Violet method (Hach Method 8028). Range: 0.002 - 0.275 mg/L LAS. Reagent Set includes: Benzene, Buffer Solution and Detergent Reagent Powder Pillows for 70 - 100 tests. Note: No Air Shipments - see Product Number 2446801 for Air Freight Configuration.

---

### Specifications

Instrument:	All except PC II
Method:	8028
Method Name:	Crystal Violet
Number of tests:	70 - 100
Range:	0.002 - 0.275 mg/L LAS

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)