



Be Right™



Phenol Reagent Powder Pillows, pk/100

Product #: 87299
SGD Price: Contact Hach
Available

For determination of phenols by the 4-Aminoantipyrine Method. Hach method 8047. Pack of 100 powder pillows.

Specifications

Description :	Phenol Reagent Powder Pillows
Method Name:	4-Aminoantipyrine
Number of tests:	50
Parameter:	Phenols
Platform :	Powder Pillows
Quantity:	100 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)