



Be Right™



Reagent Set, Glycol Test Kit

Product #: 2446900
SGD Price: Contact Hach
No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of glycol by the Hach EG-1 Glycol Test Kit (2186400). Reagent set includes: Glycol Test Reagent 1, Glycol Test Reagent 2, Glycol Test Reagent 3 and Ethylene Glycol - 0.05% solution. Approximately 100 tests.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)