



Total Organic Carbon (TOC) Reagent Set, LR

Product #: 2760345
SGD Price: Contact Hach
Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Total Organic Carbon (TOC) by the Low Range Test 'N Tube™ method.

Hach Method 10129. Range: 0.3 - 20.0 mg/L C

Set includes digestion vials and reagents for 25-50 tests (dependent upon the number of reagent blank tests performed).

Requires digestion.

Affordable

The Hach TOC method requires only a fraction of the investment and operational costs of a conventional TOC analyzer.

Easy to use

Only minimal training is required to perform the method.

Accurate

The Hach TOC method correlates well to regulatory-approved instrumental methods on natural waters, drinking water, and wastewater samples.

Non-hazardous

Prepared reagents are contained in screw-capped vials to reduce laboratory waste and associated disposal problems.

Specifications

Digestion Required:	Yes
Method:	10129
Method Name:	Direct
Number of tests:	25 - 50
Parameter:	Total Organic Carbon
Platform :	Test 'N Tube
Range:	0.3 - 20.0 mg/L C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)