



Be Right™



Carbon Dioxide Reagent Set for Model CA-23 Test Kit

Product #:	2437500
SGD Price:	Contact Hach
Available	

For determination of Carbon Dioxide by drop-count titration. Replacement reagent set for Hach Test Kit Model CA-23. Includes: Phenolphthalein Indicator Solution and Sodium Hydroxide Standard Solution (0.01 N).

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)