



Be Right™



Stablcal Stabilised Formazine Turbidity Standard #0.1 NTU

Product #: 2659749
SGD Price: Contact Hach
Available

Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, <0.1 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Specifications

Concentration:	< 0.1 NTU
Platform :	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)