



Be Right™



## Stablcal stabilised Formazine Turbidity Standard, 100 NTU, 100 mL

**Product #:** 2660242  
**SGD Price:** Contact Hach  
**Available**

Ready to use, stabilised Formazine Primary standard solution in 100 mL bottle, 100 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

### Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### Based on Formazin Primary Standards

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Specifications

|                |          |
|----------------|----------|
| Concentration: | 100 NTU  |
| Platform :     | Solution |
| Quantity:      | 100 mL   |

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)