



Be Right™



Secondary Turbidity Standards Kit for 2100 Portable Turbidimeters

Product #: 2464105

SGD Price: Contact Hach
Available

Kit includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU.

Gelex Secondary Standards Kit, for 2100P and 2100Q Portable Turbidimeters (includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU).

Specifications

Quantity: set of 3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)