



Be Right™



## SpecCheck Secondary Gel Standards Set, DPD Chlorine - HR

Product #: 2893300

SGD Price: Contact Hach  
Available

DPD Chlorine, High Range, 0-6.5 mg/L as Cl<sub>2</sub>. Used to confirm consistent instrument response. Set contains four vials filled with gels that simulate the test color at various concentrations.

\*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

### Fast and Convenient

Check instrument response without sample or reagent preparation.

### Stable

Colored gels are stable - assuring consistent performance

---

## Specifications

|                |   |
|----------------|---|
| Concentration: | 0 - 8.0 mg/L Cl <sub>2</sub>                          |
| Description :  | SpecCheck Gel Secondary Standard Kit-HR Chlorine, DPD |
| Instrument:    | Includes values for all instruments                   |
| Method Name:   | DPD   |
| Quantity:      | Set of 4 Vials  |