



Be Right™



Buffer Powder Pillows, pH 10.01 Colour-coded Blue, pk/50

Product #: 2227166
SGD Price: Contact Hach
Available

pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder pillow prepares 50 mL of fresh buffer solution. Colour-coded for easy identification (Blue). pH is 10.01 ±0.02 at 25°C. Pack of 50 Powder Pillows

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

| | |
|---------------------|-------------------------------------|
| Color coded?: | Yes |
| Description : | pH Buffer Calibration reagent |
| NIST Traceable: | Yes |
| Package Type: | Powder Pillows |
| pH Value(s): | 10.01 ±0.02 pH at 25°C |
| Quantity: | 50 pcs |
| Shelf Life: | 24 months |
| Special Remarks: | Discard prepared solution after use |
| Storage conditions: | 10 - 25 °C (protect from light) |