



Be Right™



## Buffer Solution, pH 7.00, Colour-coded Yellow, 4L

Product #: 2283556  
 SGD Price: Contact Hach  
 Available

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colour-coded pH Buffer solution for easy identification. pH is 7.00 ±0.02 at 25°C. 4 L bottle

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Color coded?:	Yes
Description :	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
pH Value(s):	7.00 ±0.02 pH at 25°C
Quantity:	4 L
Region:	EU
Shelf Life:	24 months
Special Remarks:	Observed shelf-life after opening, following recommended storage and use: 6 months
Storage conditions:	10 - 25 °C (protect from light)