



Be Right™



Buffer Solution, pH 10.01, Colour-coded Blue, 500 mL

Product #: 2283649
SGD Price: Contact Hach
Available

pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colour-coded pH Buffer solution for easy identification. pH is 10.01 ±0.02 at 25°C. 500 mL bottle

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Color coded?:	Yes
Description :	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
Parameter:	pH
pH Value(s):	10.01 ±0.02 pH at 25°C
Quantity:	500 mL
Shelf Life:	12 months
Storage conditions:	10 - 25 °C (protect from light)