



NitraVer® 5 Nitrate Reagent Powder Pillows, 25 mL, pk/1000

Product #: 1403428
SGD Price: Contact Hach
No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For mid-range and high range nitrate determination by the NitraVer® 5 cadmium reduction method. Hach methods 8039 and 8171. Range: up to 30 mg/L as N. Sample Size: 25 mL. Pack of 1000 powder pillows.

Specifications

Description :	NitraVer 5 Reagent Powder Pillows, 25 mL
Instrument:	DR/3000, DR/2000
Method:	8039
Method Name:	Cadmium Reduction
Method Number 2:	8171
Number of tests:	1000
Platform :	Powder Pillows
Quantity:	1000 /pk
Range:	0.3 - 30.0 NO ₃ ⁻ -N
Range 2:	0.1 - 10.0 mg/L NO ₃ ⁻ -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)