



Be Right™



## NitriVer® 2 Nitrite Reagent Powder Pillows, 25 mL, pk/100

**Product #:** 221969  
**SGD Price:** Contact Hach  
**Available**

For high range nitrite determination by the ferrous sulfate method. Hach method 8153. Range::2 - 150 mg/L NO<sub>2</sub>-N. Sample Size: 25 mL. Pack of 100 powder pillows.

---

### Specifications

Description :	NitriVer 2 Nitrite Reagent Powder Pillows, 25 mL
Instrument:	DR/2000, DR/3000
Method:	8153
Method Name:	Ferrous Sulfate
Number of tests:	100
Platform :	Powder Pillows
Quantity:	100 /pk
Range:	2 - 150 mg/L NO <sub>2</sub> <sup>-</sup>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)