



Be Right™



Nitrogen-Nitrate Reagent Set, LR, 25 mL, Cadmium Reduction

Product #: 2245100
SGD Price: Contact Hach
Ships within 3 days

For determination of Nitrate Nitrogen by the Cadmium Reduction method. Hach Method 8192. Range: up to 0.50 mg/L NO₃-N. Sample Size: 25 mL. Reagent Set includes: NitraVer® 3 Reagent Powder Pillows and NitriVer® 6 Reagent Powder Pillows for 100 tests.

Specifications

Description :	Nitrate Reagent Set, 25 mL, LR
Instrument:	DR/3000, DR/2000
Method:	8192
Method Name:	Cadmium Reduction
Number of tests:	100
Platform :	Powder Pillows
Range:	0 - 0.50 mg/L NO ₃ ⁻ -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)