



Be Right™



Nitrogen-Ammonia Reagent Set, Nessler

Product #: 2458200
SGD Price: Contact Hach
Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

For determination of Ammonia Nitrogen by the Nessler Method. USEPA accepted for reporting wastewater analysis (distillation required) using Hach Method 8038. Range: 0.02-2.50 mg/L NH₃-N. Reagent set contains Nessler Reagent, Mineral Stabilizer and Polyvinyl Alcohol Dispersing Agent for approximately 250 tests.

Specifications

EPA compliant:	Yes
Instrument:	All except DR 900, DR/800 Series, PC II
Method:	8038
Method Name:	Nessler
Number of tests:	250
Parameter:	Ammonia
Platform :	Solution
Range:	0.02 - 2.50 mg/L NH ₃ -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)