



Be Right™



NitraVer® 6 Nitrate Reagent Powder Pillows, 30 mL, pk/1000

Product #: 1411928

SGD Price: Contact Hach
No availability shown

For determination of low range nitrate by the cadmium reduction method. Range: up to 0.50 mg/L as nitrate-N. Use with NitriVer® 3. Sample Size: 30 mL.

Specifications

Method Name: Cadmium Reduction
Platform : Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)