



Be Right™



NitriVer 3 Nitrite Reagent Powder Pillows, 5 mL, pk/100

Product #: 1407899
SGD Price: Contact Hach
Available

For nitrite determination by the NitriVer 3 diazotization method. Sample size 5 mL. For use with Hach visual test kits. Range: up to 0.5 mg/L NO₂-N.

Specifications

Method Name: Diazotization
Platform : Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)