



Be Right™



NitraVer® 6 Nitrate Reagent Powder Pillows, 5 mL, pk/100

Product #:	1412099
SGD Price:	Contact Hach
Available	

For determination of low range nitrate by the cadmium reduction method. For use with Hach visual test kits. Use with NitriVer® 3. Sample Size: 5 mL.

Specifications

Method Name:	Cadmium Reduction
Platform :	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)