



Be Right™



Nitrogen-Nitrate Reagent Set, LR, 5 mL Sample

Product #: 2460800
SGD Price: Contact Hach
Available:

For determination of Nitrogen-Nitrate by the Cadmium Reduction method. Sample Size: 5 mL For visual test kits. Reagent set includes: NitriVer® 3 Powder Pillows and NitraVer® 6 Reagent Powder Pillows for 100 tests.

Specifications

Method Name: Cadmium Reduction

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)