



Be Right™



Boric Acid Solution (Indicating), 20 g/L, 1 L

Product #: 2339053
SGD Price: Contact Hach
Available

For determination of ammonia by the phenate method. 20 g/L H₃BO₃, with indicators, APHA formulation 1 liter bottle

Specifications

Platform : Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)