



Be Right™



## Aqua Ammonia Reagent Set, Digital Titrator

**Product #:** 2930500  
**SGD Price:** Contact Hach  
**No availability shown**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Set for Model NI-HRDT Digital Titrator Kit, Cat. No. 2930400. Range: 5-35% (50-350 g/L) NH<sub>3</sub>. Includes Sulfuric Acid Digital Titrator Cartridge and Wide Range pH Indicator. Approximately 100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifications

Method:	10222
Number of tests:	100
Platform :	Digital Titrator
Range:	5 - 35 % = 50-350 g/L NH <sub>3</sub>

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)