



Be Right™



Ammonia Standard Solution, 10 mg/L, 500 mL

Product #: 15349
SGD Price: Contact Hach
Available

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Ammonia, 10.0 mg/L \pm 0.5 mg/L NH₃-N. 500 mL bottle

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	10.0 mg/L \pm 0.5 mg/L NH ₃ -N
Concentration 2:	0.00435 mol/L \pm 0.00003 mol/L Na ⁺
Package Type:	Bottle
Parameter:	Ammonium
Platform :	Solution
Quantity:	500 mL
Shelf Life:	5 years
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)