



Be Right™



NitriVer® 3 Nitrite Reagent Powder Pillows, 10 mL, pk/100

Product #: 2107169
SGD Price: Contact Hach
Available

For determination of low range nitrite by the NitriVer® 3 diazotization method.

USEPA accepted for reporting wastewater analysis using Hach Method 8507.

Range: up to 0.30 mg/L NO₂-N. Sample Size:10 mL.

Pack of 100 Powder Pillows

Specifications

Description :	NitriVer 3 Nitrite Reagent, 10 mL
EPA compliant:	Yes
Instrument:	All except PC II
Method:	8507
Method Name:	Diazotization
Number of tests:	100
Parameter:	Nitrite
Platform :	Powder Pillows
Quantity:	100 /pk
Range:	0.002 - 0.300 mg/L NO ₂ ⁻ -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)