



Be Right™



NitraVer® 6 Nitrate Reagent Powder Pillows, 10 mL, pk/100

Product #: 2107249
SGD Price: Contact Hach
Available

For determination of low range nitrate by the NitraVer® 6 cadmium reduction method. Range: 0.01 to 0.50 mg/L as nitrate-N. Use with NitriVer® 3. Sample Size: 10 mL. Pack of 100 Powder Pillows

Specifications

Description :	NitraVer 6 Powder Pillows
Method Name:	Cadmium Reduction
Number of tests:	100
Parameter:	Nitrate
Platform :	Powder Pillows
Quantity:	100 /pk
Range:	0.01 - 0.50 mg/L NO_3^- -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)