



Be Right™



## Ammonia Salicylate Reagent Powder Pillows, 10 mL, pk/100

**Product #:** 2653299  
**SGD Price:** Contact Hach  
**Available**

For determination of Ammonia Nitrogen by the Salicylate method. Hach Method 8155.

Range: 0.01 - 0.50 mg/L NH<sub>3</sub>-N. Sample Size: 10 mL.

Pack of 100 Powder Pillows.

---

### Specifications

Method:	8155
Method Name:	Salicylate
Number of tests:	100
Parameter:	Ammonia
Platform :	Powder Pillows
Quantity:	100 /pk
Range:	0.01 - 0.50 mg/L NH <sub>3</sub> -N

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)