



Be Right™



Nitrate Accuvac, 0.3-30.0 or 0.1-10.0 mg/L NO₃-N

Product #: 2511025
SGD Price: Contact Hach
Available

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Nitrogen, Nitrate by the cadmium reduction method. Pack of 25.

Fast
AccuVac® Ampules provide quick results. Just snap, fill and read - saving you valuable time.

Easy
No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

Accurate
Optical quality evacuated vials assure precise measurements and reagent stability.

Use them anywhere you need results.
Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

Specifications

Instrument:	All
Method:	8039
Method Name:	Cadmium Reduction
Method Number 2:	8171
Number of tests:	25
Parameter:	Nitrate
Platform :	Accuvac
Range:	0.3 - 30.0 mg/L NO ₃ -N
Range 2:	0.1 - 10.0 mg/L NO ₃ -N

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)