



Be Right™



NitraVer X Nitrogen-Nitrate Reagent Set, HR

Product #: 2605345
SGD Price: Contact Hach
Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of High Range Nitrate Nitrogen by the Chromotropic Acid Test 'N Tube Method. Hach Method 10020. Range: 0.2 - 30.0 mg/L NO₃-N

50 tests per package.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Specifications

Description :	NitraVer X Nitrate Reagent Set
Method:	10020
Method Name:	Chromotropic Acid
Number of tests:	50
Parameter:	Nitrate
Platform :	Test 'N Tube
Quantity:	50 /pk
Range:	0.2 - 30.0 mg/L NO ₃ -N
Special Remarks:	Set contains: Cuvette test tubes, NitraVer X Reagent A (2605400), NitraVer X Reagent B Powder Pillows (2605546)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)