



Total Nitrogen Reagent Set, HR, TNT

Product #: 2714100
SGD Price: Contact Hach
Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Total Nitrogen by the Persulfate Digestion Test 'N Tube method:

Hach Method 10072. Range: 10 - 150 mg/L N.

25-50 tests (dependent upon the number of reagent blank tests performed)

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Perform digestions and photometric measurements in the same vial.

Save time, space and effort.

Specifications

Description :	Nitrogen, Total, TNT Reagent Set, HR
Digestion Required:	Yes
Footnote:	†Dependent on the number of performed reagent blanks.
Method:	10072
Method Name:	Persulfate Digestion
Number of tests:	25 - 50
Parameter:	Nitrogen, total
Platform :	Test 'N Tube
Range:	2 - 150 mg/L N
Special Remarks:	Set consists of 2 cartons: 27140-45 and 26721-45 27140-45 = TN Hydroxyde reagent solution (2714045), Persulfate Powder Pillows (2671846) 26721-45 = TN Acid reagent C solution (2672145), Bisulfite reagent A Powder Pillows (2671946), Indicator reagent B Powder Pillows (2672046)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)