



**Be Right™**



## **MultiTest kit case (218 x 246 x 437 mm), blue polypropylene**

**Product #:** 4660700  
**SGD Price:** Contact Hach  
**Ships within 1 week**

Durable blue polypropylene carrying case allows easy transport of items for field testing.

**Recessed handle in lid provides easy stacking**

**Contents are carried in upright position**

**Includes recessed flat-panel areas for labeling**

---

### **Specifications**

Case Style: Q  
Dimensions (H x W x D): 218 mm x 246 mm x 437 mm  
Height : 8.5 in