



**Be Right™**



## **Large custom kit case assembly (574 x 325 x 295 mm)**

<b>Product #:</b>	<b>4766600</b>
SGD Price:	Contact Hach
Available	

Designed to withstand the harshest conditions, this rugged, blue polypropylene case assembly includes a foam insert to prevent kit components from shifting, breaking or spilling.

**Allows easy transport of instruments, reagents, and apparatus into the field**

**Chemical-resistant test kit case is molded from durable, easy-to-clean polypropylene**

---

### **Specifications**

Dimensions (H x W x D): 574 mm x 325 mm x 295 mm