



Be Right™



Polypropylene Low Form Beaker, 600 mL

Product #: 108052
SGD Price: Contact Hach
Available

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Non-toxic, noncontaminating, and biologically inert

Double spout for easy pouring

Autoclavable

Specifications

Capacity:	600 mL
Diameter:	110.7 mm with spout
Instrument:	All except DR/4000, PC II
Method:	8006
Quantity:	1 each