



Be Right™



pH Test Kit, 6.6 - 8.4 pH, Model 17H

Product #: 147008
SGD Price: Contact Hach
Available

Simple, low-cost color disc method is the most accurate visual test kit method.

Achieve fast results with an easy-to-carry test kit

Achieve fast results with an easy-to-carry test kit

Specifications

| | |
|---------------------|-----------------------|
| Case Style: | B |
| Method Name: | Color disc/Phenol red |
| Model: | 17H |
| Number of tests: | 200 |
| Parameter: | pH, mid range |
| Platform : | Color Disc |
| Power supply: | NA |
| Range: | 6.6 - 8.4 pH |
| Range 2: | NA |
| Range 3: | NA |
| Ship Wt. (lbs): | 1 |
| Smallest Increment: | 0.1 pH |

What's in the box?

Includes indicator reagent, color disc, comparator box, two viewing tubes, instructions and carrying case.