



Be Right™



Hardness (Total) Test Kit, Model HA-DT

Product #: 2063600
SGD Price: Contact Hach
Available

Utilizing the Hach Digital Titrator, the Model HA-DT kit is used to determine hardness concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifications

Case Style:	F
Method Name:	Digital Titrator/EDTA
Model:	HA-DT
Number of tests:	100
Parameter:	Hardness, total - As CaCO_3
Platform :	Digital Titrator
Range:	10 - 4,000 CaCO_3
Ship Wt. (lbs):	4
Smallest Increment Steps:	0.1-10

What's in the box?

Includes four reagents, flask, digital titrator, delivery tube, cylinder, instruction sheet, and carrying case.