



Be Right™



Chlorine Dioxide Color Disc Test Kit, Model CLD-2

Product #: 2842800
SGD Price: Contact Hach
Available

Monitor ClO₂ levels on-site with the Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Specifications

Case Style:	D
Method Name:	Color Disc/DPD
Model:	CLD-2
Number of tests:	100
Parameter:	Chlorine Dioxide
Platform :	Color Disc
Power supply:	NA
Range:	0 - 2.0 ClO ₂
Range 2:	NA
Range 3:	NA
Ship Wt. (lbs):	3
Smallest Increment:	0.1

What's in the box?

Includes two reagents, tube, color comparator box, color discs, instruction sheets, and carrying case.