



Total Chlorine Test Kit, Model CN-DT, 20-2000 mg/L

Product #: 2471100
SGD Price: Contact Hach
Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Kit contains all required reagents and apparatus to perform analysis in the field

Specifications

Case Style:	F
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Method Name:	Digital Titrator/Thiosulfate
Model:	CN-DT
Number of tests:	100
Parameter:	Chlorine, mid and high range - As total Cl_2
Platform :	Digital Titrator
Range:	20 - 2,000 mg/L Cl_2
Ship Wt. (lbs):	4
Smallest Increment:	0.2 - 5.0
Smallest Increment Steps:	0.2 5.0

What's in the box?

Includes four reagents, flask, digital titrator, delivery tube, cylinder, bottle, clippers, instruction sheet, and carrying case.
Product details pdf footer