



Be Right™



## Dissolved Oxygen AccuVac® Kit

**Product #:** 2515050  
**SGD Price:** Contact Hach  
**Ships within 3 days**

Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifications

Case Style:	G
Method Name:	Color disc/HRDO
Model:	Color disc/HRDO
Number of tests:	25
Parameter:	Dissolved Oxygen - As DO
Platform :	Accuvac
Power supply:	NA
Range:	0 - 15 mg/L O <sub>2</sub>
Range 2:	NA
Range 3:	NA
Ship Wt. (lbs):	4.1
Smallest Increment:	1

---

### What's in the box?

Includes one reagent, beaker, color comparator box, vial, beaker, ampule breaker, color discs, instruction sheet, and carrying case.