



Be Right™



## Alkalinity Test Kit, Model AL-AP, gpg

**Product #:** 2444300  
**SGD Price:** Contact Hach  
**Available**

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.  **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Method Name:	Drop count titration/Sulfuric acid
Model:	AL-AP
Number of tests:	100
Parameter:	Alkalinity - (P) & Total (MO) as CaCO <sub>3</sub>
Platform :	Drop Count Titrator
Range:	- 8 gpg
Range 2:	1 - 20 gpg**
Ship Wt. (lbs):	1
Smallest Increment:	0.4 gpg  1 gpg

### What's in the box?

Includes sample measuring tube, mixing bottle, reagents, instructions and carrying case.