



Be Right™



Alkalinity Test Kit, Model AL-AP, mg/L

Product #: 2444301
 SGD Price: Contact Hach
 Available

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Method Name:	Drop count titration/Sulfuric acid
Model:	AL-AP, MG-L
Number of tests:	100
Parameter:	Alkalinity - (P) & Total (MO) - as CaCO ₃
Platform :	Drop Count Titrator
Range:	5 - 100 mg/L
Range 2:	20 - 400 mg/L
Ship Wt. (lbs):	1
Smallest Increment:	5
	20

What's in the box?

Includes sample measuring tube, mixing bottle, reagents, instructions and carrying case.