



Be Right™



Ozone Test Kit, Model OZ-2

Product #: 2064400
SGD Price: Contact Hach
Available

This compact and low cost ozone test kit utilizes an easy to use color disc method. Simply grab a sample, add powder and compare sample color to the color wheel to arrive on an accurate measurement.

Hach's most economical ozone test kit

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with a compact, easy-to-carry test kit

Specifications

Case Style:	B
Method Name:	Color disc/DPD
Model:	OZ-2
Number of tests:	100
Parameter:	Ozone - As O ₃
Platform :	Color Disc
Power supply:	NA
Range:	0 - 2.2 mg/L O ₃
Range 2:	NA
Range 3:	NA
Ship Wt. (lbs):	1
Smallest Increment:	0.05

What's in the box?

Includes color disc, comparator box and tubes, reagents, instructions and carrying case.