



Be Right™



## Phosphonate Test Kit, Model PN-10, without UV Lamp and Power Supply

Product #: 2113301  
 SGD Price: Contact Hach  
 Ships within 3 days

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests. Purchase UV lamp and power supply for the appropriate voltage separately.

**Simple, fast, color disc method**

**UV lamp digestion converts phosphonates into reactive phosphate**

**Control phosphonate concentration to save chemical costs**

**Hach's reliable PhosVer 3 reagent provides results in under 15 minutes**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specifications

Case Style:	F
Method Name:	PhosVer 3/Color disc
	(Does not include UV lamp or power supply.)
Model:	PN-10
Number of tests:	0 - 100
Parameter:	Phosphonate - As $PO_4^{3-}$
Platform :	Color Disc
Power supply:	115 Vac UV Lamp and Power Supply
Range:	0 - 100 mg/L $PO_4$
Range 2:	NA
Range 3:	NA
Ship Wt. (lbs):	5
Smallest Increment:	1

### What's in the box?

Includes color disc, comparator box and tubes, volumetric flask, graduated cylinder, reaction bottle, reagents, instructions and carrying case. Purchase UV lamp and power supply for the appropriate voltage separately.