



Be Right™



Immunoassay Reagent Set, Atrazine in Water Pocket Colorimeter II Analysis System, 100 Tests

Product #:

2762710

 **Hazardous**

SGD Price:

Contact Hach

Call for ship date

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Atrazine in water by the Immunoassay Method.

What's in the box?

Includes atrazine enzyme conjugate, color developing solution, stop solution and three calibrator solutions. For approximately 100 tests.