



Be Right™



Vacuum Pump, 3.5 CFM, 60 Hz

Product #:	2947000
SGD Price:	Contact Hach
Available	

"Provides up to 27.6" Hg vacuum for applications in laboratory, office or industrial plant—vacuum distillation, filtration, and air pollution studies employing membrane

filters. Quiet, pulse less air delivery and vacuum. Complete with vacuum and pressure gauges and regulator, hose barbs, intake jar assembly, muffler, switch, cord, and plug. Capacity: 128 liters/min (4.5 CFM). 1/3 HP. Dimensions (LxWxH): 15x10x10 inches. Weight: 16.5 lbs. Tubing needed: 3/8"."

Specifications

Electrical Specifications:	115 Vac, 60 Hz
Length:	381 mm
Quantity:	1 each